

NOTICE: This document contains references to Varian. Please note that Varian, Inc. is now part of Agilent Technologies. For more information, go to [www.agilent.com/chem](http://www.agilent.com/chem).



## APPLICATION NOTE

Application 1928 - HPLC

### Ginsenosides

#### Ginseng

Column: Polaris™ 5μ C18-A, 150 x 4.6mm  
Part nr.: A2000150X046  
Eluent: A: MeCN:H<sub>2</sub>O 15:85 B: MeCN:H<sub>2</sub>O 80:20  
Gradient: 0 to 30%B in 70 min.  
Flow Rate: 1.5 ml/min  
Detection: UV @ 205nm  
Sample: 1. Rg1  
2. Re  
3. Rf  
4. Rb1  
5. Rbc  
6. Rb2  
7. Rd

